


Tutorial:

How to train your own StarDist 3D model using google colab

Please first follow the tutorial how to generate 3D training samples using LabKit to train a StarDist model.

- 1) Open you google drive account: <https://drive.google.com/drive/my-drive>
- 2) Import the input images and masks folder



My Drive > Training_mydataset_02-11-2020

Name ↑	Owner	Last modified	File size
images	me	2:29 PM me	—
masks	me	2:29 PM me	—

Images must be .tiff files and input images and masks **must have the same filename**. Remember that the input images must only have the DAPI channel.

Example:

Input images:

C:\Users\fabbe\Desktop\MyDataSet\Training\images\Image00.tif

C:\Users\fabbe\Desktop\MyDataSet\Training\images\Image01.tif

Masks:

C:\Users\fabbe\Desktop\MyDataSet\Training.masks\Image00.tif

C:\Users\fabbe\Desktop\MyDataSet\Training.masks\Image01.tif

- 3) Open the google colab notebook (if the link is dead please go to the github link below):
https://colab.research.google.com/github/HenriquesLab/ZeroCostDL4Mic/blob/master/Colab_notebooks/StarDist_3D_ZeroCostDL4Mic.ipynb#scrollTo=4waLStmORPFo


Additional informations about it and acknowledgements:

<https://github.com/HenriquesLab/ZeroCostDL4Mic>


- 4) Go through the notebook as explained, some additional information about each step:

Step 3.1

To go through your google drive file click the folder icon as shown below



- 1) Navigate through your folders to copy the path containing the Training_source (images of nuclei) and Training_target (masks).
- 2) The model path is the path where the model folder will be created.
- 3) The model name is the folder name that will be created for your model, you don't need to create manually this folder it will be created automatically.
- 4) Don't forget to uncheck the Use_Default_Advanced_Parameters cell!
- 5) Set the percentage_validation as 0
- 6) Recommendation: set n_rays as 192
- 7) Set the patch_size (xy) and patch_heigh (z) as big as possible (it shouldn't be larger and deepest than your smallest image)



Step 3.2

You should use the data augmentation, don't forget to activate the Use_Data_augmentation cell.




- 1) Deactivate the elastic transformation cell
- 2) Activate the rotation cell
- 3) Activate the flip cell
- 4) Deactive the save_augmented_images cell

Step 3.3

If you train your model from scratch, deactivate the Use_pretrained_model cell. Otherwise if you have already trained a model and want to refine it, considers training from your pretrained model as it can lead to better results.





If you decide to train a model from a previous model, put your model in google drive as show below and copy the path in the pretrained_model_path cell.

My Drive > Training_mydataset_02-11-2020 ▾

Name ↑	Owner	Last modified	File size
 images	me	2:29 PM me	—
 masks	me	2:29 PM me	—
 <u>old_model</u>	me	Aug 13, 2020 me	—

The folder model should contain the following files:

My Drive > Training_mydataset_02-11-2020 > old_model ▾

Name ↓	Owner	Last modified	File size
 weights_now.h5	me	Aug 4, 2020 me	6 MB
 weights_best.h5	me	Aug 4, 2020 me	6 MB
 thresholds.json	me	Aug 4, 2020 me	26 bytes
 config.json	me	Aug 4, 2020 me	1 KB

Step 4.2

Once done you can download your new model in your google drive in the folder model name located in model path (see steps 3.1)